

## Experiment No. 7

### Determination of $K_La$

#### Objective:

To determine  $K_La$  in water

#### Introduction:

In the model system used for the experiment to determine  $K_La$  in water i.e., in the case where no reaction is taking place, the following simplifications are valid:

- $R_0 = 0$  because there is no  $O_2$  sink in the system
- Only data for uptake will be considered in the following equation:

$$\frac{dC_0}{dt} = OTR - OUR = K_La(C_0^* - C_0) - R_0$$

The concentration  $C_0$  can be achieved by flushing nitrogen in the system. At time  $t_0$  the degassing with nitrogen is stopped and from time  $t_0$  there is constant aeration. The next phase is one of uptake. The curve obtained for  $C_0 = f(t)$  is followed until there are sufficient constant values at the upper end of the curve. For this event following equation is valid:

$$\frac{dC_0}{dt} = OTR = K_La(C_0^* - C_0)$$

Solving the above differential equation for the starting conditions  $C_0(t_0 - t) = C_0^0$  results in equation:

$$C_0(t) = C_0^* - (C_0^* - C_0^0) \exp[-K_La(t - t_0)]$$

To determine the  $K_La$  value by the dynamic method, short bursts of measurement with an electrode are necessary so that measurement time with the electrode has no influence on the value of  $K_La$ .

#### Media and materials:

Tap water is used as the medium. Aeration is provided by an aeration system below the stirrer. Compressed air or nitrogen to be used for aeration or degassing the system. The fermenter is thermostatically controlled.

**Procedure:**

1. Calibrate the DO probe according to the manufacturers' instruction
2. Set the temperature of the fermenter to 20°C
3. Position the DO probe so that the membrane is at the height of the stirrer. Set the probe to measure the percent DO level and allow the probe to reach the temperature of the fermenter.
4. Switch on and set data collector to measure every 5 seconds.
5. Perform the experiment with 1.5 l of medium (water) at 250 rpm
6. Repeat the experiment at various agitation rates of 500, 750 and 1,000 rpm.

**Task Required**

**Calculation of  $K_La$  as a function of agitator speed.**